



Easy to clean advanced polymer material is a marked improvement over traditional metal carts: won't dent, chip, rust, flake or corrode.

Molded ergonomic handle is easy to grip.

Interchangeable 3", 6" and 9" (76, 152 and 230mm) drawers with removable totes can be fully extended for easy access.

Side storage allows easy customization of a wide variety of accessories from trash can and Sharps container to storage bins and shelves.

The convenient overbridge can be fitted with a wide selection of baskets, shelves, and bins to keep necessary items within reach.

Swing-out Side Storage units increase work surface up to 135%.

Color-coded drawer pulls are available in an array of choices to fit any system or decor.

Label holders available to quickly identify drawer contents.

5" (127mm) polyurethane casters with easily accessed toe brake make for easy rolling.



Starsys Carts

Starsys provides a complete system of enclosed carts in a wide variety of heights and system widths. When mobility is a key part of your storage needs, Starsys carts provide a broad selection of unique solutions. This modular system approach allows you to create a cart to meet your exact needs.

Starsys cart widths:

The Starsys cart system is built upon a square module. Multiplying this storage module results in system widths of single 21.5" (546mm), double 42" (1067mm) and triple widths 63" (1600mm).

Starsys cart heights:

Starsys carts are available in the following working heights: 34" (864mm), 37" (940mm), 40" (1016mm), 43" (1092mm) and 46" (1168mm).

General Overview:

Starsys carts feature a 5" (127mm) swivel caster with a polymer wheel as a standard. The front two casters on a Starsys cart are equipped with color-matched toe brakes.

Starsys carts feature a smooth polymer top/worksurface. If chemical resistance is a concern or heavy equipment is likely to be placed on the top of the cart, Starsys Mobile WorkCenters may be a more appropriate product selection.

An extended mobile base is required when Heavy-Duty drawers, Active Level shelving or qwikSLOT shelving is specified in a Starsys configuration. The extended mobile base adds an additional 3" (76mm) to the overall depth of the cart. The extended base combines additional counter weights with a deeper footprint to reduce any chance of overbalance when active level shelving is fully extended.

To configure a Starsys unit to your exact specifications, please visit www.metroconfigurator.com

STARSYS™ PROCEDURE CARTS



Complete configuration details are available at:
www.metro.com/catalog. Click on "online catalog."

40"

Close-up
of system



General Supply
 Cat. No. **SXS40CM1**
 List Price 2,445.00 each



Isolation
 Cat. No. **SXSISO**
 List Price 2,653.50 each



Med Surg
 Cat. No. **SXSMDSRG**
 List Price 2,586.00 each



Dressing
 Cat. No. **SXSDRS**
 List Price 2,739.00 each



Laptop Ready Cart
 (Keyboard not included)
 Cat. No. **SXSK1120L**
 Call for pricing



Laptop Ready Cart
 (Keyboard not included)
 Cat. No. **SXSK2210L**
 Call for pricing



Laptop Ready Cart
 (Keyboard not included)
 Cat. No. **SXSK4110L**
 Call for pricing

43"



Bedside Cart
 Cat. No. **SXS43CM1**
 List Price 2,304.00 each



Nurse Server
 Cat. No. **SXS43CM3**
 List Price 2,564.00 each



General Supply
 Cat. No. **SXS43CM4**
 List Price 3,831.00 each



General Supply
 Cat. No. **SXS43CM5**
 List Price 3,034.00 each



General Supply
 Cat. No. **SXS43CM6**
 List Price 3,895.00 each



Difficult Airway
 Cat. No. **SXS AIR1**
 List Price 4,597.00 each



Anesthesia
 Cat. No. **SXSANES**
 List Price 4,524.50 each



Bedside
 Cat. No. **SXS BED**
 List Price 3,897.50 each

Starsys Carts are built to exact customer specifications, therefore, they are non-cancelable, non-returnable.

Complete configuration details are available at: www.metro.com/catalog. Click on "online catalog."

43"



Decentralized Medication
Cat. No. **SXSDCMED**
List Price 4,682.50 each



Arthroscopic
Cat. No. **SXSARTH**
List Price 4,087.50 each



General Procedure
Cat. No. **SXS43CM2**
List Price 4,333.00 each



I.V. Therapy
Cat. No. **SXSIV**
List Price 4,375.00 each



Laparoscopic
Cat. No. **SXSLAP**
List Price 4,272.00 each



Unit Dose Medication
Cat. No. **SXSUDOSE**
List Price 5,463.00 each



Trauma
Cat. No. **SXSTRAMA**
List Price 4,571.50 each

46"



Biomedical
Cat. No. **SXS BIO**
List Price 4,542.50 each



Critical Care
Cat. No. **SXS CCU**
List Price 4,237.50 each



Cast
Cat. No. **SXS CAST**
List Price 4,061.00 each

To configure a Starsys unit to your exact specifications, please visit www.metroconfigurator.com



Complete configuration details are available at:
www.metro.com/catalog. Click on "online catalog."

Basic Single Wide, Locking Carts

(in.)	(mm)	Drawer Configuration	Cat. No.	List Price Each
24 ⁷ / ₈ x22 ¹¹ / ₁₆ x39 ³ / ₄	632x576x1010	1-3", 3-6", 1-9", 0-12"	SXS1310L	2,526.00
24 ⁷ / ₈ x22 ¹¹ / ₁₆ x39 ³ / ₄	632x576x1010	3-3", 2-6", 1-9", 0-12"	SXS3210L	2,616.40
24 ⁷ / ₈ x22 ¹¹ / ₁₆ x39 ³ / ₄	632x576x1010	5-3", 1-6", 1-9", 0-12"	SXS5110L	2,706.80
24 ⁷ / ₈ x22 ¹¹ / ₁₆ x42 ³ / ₄	632x576x1086	1-3", 2-6", 2-9", 0-12"	SXS1220L	2,600.00
24 ⁷ / ₈ x22 ¹¹ / ₁₆ x42 ³ / ₄	632x576x1086	4-3", 2-6", 1-9", 0-12"	SXS4210L	2,777.80
24 ⁷ / ₈ x22 ¹¹ / ₁₆ x42 ³ / ₄	632x576x1086	6-3", 1-6", 1-9", 0-12"	SXS6110L	2,868.20
24 ⁷ / ₈ x22 ¹¹ / ₁₆ x45 ³ / ₄	632x576x1162	0-3", 3-6", 2-9", 0-12"	SXS0320L	2,724.00
24 ⁷ / ₈ x22 ¹¹ / ₁₆ x45 ³ / ₄	632x576x1162	2-3", 2-6", 2-9", 0-12"	SXS2220L	2,814.00
24 ⁷ / ₈ x22 ¹¹ / ₁₆ x45 ³ / ₄	632x576x1162	3-3", 3-6", 1-9", 0-12"	SXS3310L	2,901.95



SXS3210L



SXS3310L

Carts shown in chart have taupe pulls. Images are shown with blue pulls.

Drawer Pulls (one required per drawer)

To specify color for drawer pull, order following catalog numbers for each drawer:



Dark Taupe
SXDP-TP
List Price
3.90 each



Slate Blue
SXDP-BL
List Price
3.90 each



Jade
SXDP-GR
List Price
3.90 each



Mauve
SXDP-RD
List Price
3.90 each



Yellow
SXDP-YL
List Price
3.90 each



Code Blue
SXDP-CB
List Price
3.90 each



Black
SXDP-BK
List Price
3.90 each



Red
SXDP-RE
List Price
3.90 each



Electronic TouchPad

Security Options

Electronic TouchPad (Keyless Drawer Lock)

For use with low profile, single wide, locking carts.
Three user-programmable PIN numbers.
Programmable automatic relock function.
Programmable access-denied function.
Cat. No. **SXSTP-B** List Price 1,198.00 each

Contact your Metro representative for more information.



Electronic Push Button Lock

Electronic Push Button Locks

The electronic push button lock is a simple battery operated lock.
It has a ten-digit keypad with ten user programmable PIN numbers.

Note: Electronic Push Button Lock takes up the top 3" of cart space and comes standard with key override.



Tamper-Evident Locks

Tamper-Evident Locks

Description	Cat. No.	List Price Each
24" (610mm) Security Bar, Drawer Tabs & Seals	SXPSLK-24	149.50
27" (686mm) Security Bar, Drawer Tabs & Seals	SXPSLK-27	149.50
30" (762mm) Security Bar, Drawer Tabs & Seals	SXPSLK-30	149.50
33" (839mm) Security Bar, Drawer Tabs & Seals	SXPSLK-33	149.50
36" (914mm) Security Bar, Drawer Tabs & Seals	SXPSLK-36	149.50

Starsys Carts are built to exact customer specifications, therefore, they are non-cancelable, non-returnable.